



## Does HSQ/PQ software need to be saved as a macro enabled XLSM file - FAQ

YES, the HSQ/PQ software must be macro enabled, the file extension is XLSM.

One small downside of using Excel, lies in its reliance on macro-enabled spreadsheets. Some companies impose stringent restrictions on the use of such macros.

If your laboratory SOP requires the ability to run a spreadsheet, you may need to consult with your IT department to enable this functionality.

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 07-22-2024