



## Does the Mobile Phase Filter made with polyethylene contain PTFE - FAQ

### Do the polyethylene Mobile Phase Filters contain PTFE?

These products do not contain PTFE as long as you use a Tefzel / [ETFE ferrule](#), catalog number [58847-01](#). [Mobile Phase Filter, 20 um, Polyethylene 58847-01](#)

*[Replacement, Mobile Phase Filter, 20 um, Polyethylene. 58847-02](#)*

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 07-22-2024