



Does your virgin PEEK tubing have additives of any kind – FAQ

Our PEEK tubing does not have additives, but you may want to consider that PEEK is produced using phosphoric acid. It is possible that some remnants of this acid may be present in our PEEK Tubing.

NOTE: PEEK is not a fluoropolymer and is not impervious to all substances. It's crucial to refer to the list of materials that can potentially harm PEEK to ensure that you're not employing a solvent that could degrade it.

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 10-22-2024