

## FlexiTube specifications – Tech Information

**The FlexiTube™ length and coil data for selecting the right tube for you.**

Ordering information for each type of FlexiTube can be found on the [ordering information page](#), under each individual item.

Type	Coil resting position	Coil extended position	Coil quantity	FlexiTube length fully extended	Total length of tubing
2	10 mm	30 mm	4	330 mm	710 mm
3	20 mm	120 mm	9	420 mm	945 mm
4	35 mm	200 mm	15	500 mm	1225 mm
5	50 mm	300 mm	23	600 mm	1585 mm
6	100 mm	400 mm	46	700 mm	2060 mm
7	150 mm	500 mm	69	900 mm	3150 mm
8	170 mm	800 mm	78	1100 mm	4710 mm
9	225 mm	1100 mm	104	1400 mm	6100 mm
10	250 mm	1300 mm	116	1600 mm	6930 mm

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 10-13-2024