



## Is PTFE tubing suitable for VHP sterilization (CIP, SIP)- FAQ

Yes, our ARE-Applied Research™ brand of PTFE tubing works well with this type of sterilization given the temperature limit of 450° F.

- **NOTE: Important facts about PTFE tubing sterilization procedures:**
- *VHP is vaporized hydrogen peroxide and is a commonly used method for achieving a high level of sterility assurance.*
- *CIP-SIP systems involve the use of pure steam to ensure that large scale processes are free from impurities and organic contaminants.*

**ARE-APPLIED RESEARCH™**

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 07-27-2024