

## Can I backflush the Cogent column filter for longer in-service life - FAQ

### The Cogent column filter can be backflushed but...

Backflushing might push the internal filter element out of place. And not having a clog in the frit, but with slightly dislodged filter element suggests that filter is not working. Therefore, backflushing is not supported.

Click [HERE](#) for ordering information and a picture of the Cogent column filters.



Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 07-22-2024