



## Is there a difference in the properties of different colors of rubber in vial cap septa - FAQ

Different colors do not always indicate a difference in the septa properties of vial cap septa.

Septa can have slightly different properties regardless of color but the differences will always be indicated in the product Description or the Additional Information on the product page.

[Click HERE for Autosampler Vial Cap Ordering Information.](#)

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 05-21-2024