

Different colors do not always indicate a difference in the septa properties of vial cap septa.

Septa can have slightly different properties regardless of color but the differences will always be indicated in the product Description or the Additional Information on the product page.

[Click HERE for Autosampler Vial Cap Ordering Information.](#)

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com