



Standard silicone pre slit sealing mat compared to solvent resistant mat - FAQ

Comparison of 96 well sealing mats, silicone rubber v. chemical resistant U-2D mats.

Silicone rubber mats should only be used with pure aqueous and buffer solutions.

The solvent resistant mat can be used with methanol, acetonitrile and water / buffers. These should not be used with non-polar solvents or aggressive acid or base such as ethyl acetate.

[U-2D Silicone Mats Product Page](#)

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com