


Valganciclovir Hydrochloride USP assay with HPLC – AppNote


Alternate column used for a USP assay method of Valganciclovir Hydrochloride

This similar column for use in analyzing USP Valganciclovir Hydrochloride is easily performed, and the 10 injections overlay below demonstrates run to run consistency and great peak symmetry. The peak efficiency is very good and peak tailing guidelines are easily met with RSD values of less than 1.0%.

This demonstrates a good alternate column for your Valganciclovir Hydrochloride USP method.



Peak:



Diastereomers Valganciclovir Hydrochloride

Method Conditions:

Column: Cogent RP C18™, 5µm, 100Å

Catalog No.: 68518-15P

Dimensions: 4.6mm x 150mm

Note 1: For the full method conditions and details, consult the official United States Pharmacopeia-National Formulary (USP-NF.)

Note 2: Valganciclovir Hydrochloride is an antiviral prescription medicine for the treatment of cytomegalovirus

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

MICROSOLV

retinitis (CMV retinitis) in adults with AIDS. It is also used for the prevention of CMV disease in recipients of organ transplants who are at risk for CMV diseases.



Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com