



## What can I use for a reference solution against Holmium Oxide - FAQ

Water is an appropriate reference solution. Since the lowest wavelength in Wavelength Calibration Solution is 241 nm, any substance that doesn't absorb in that region can be used as a reference solution.

Click [HERE](#) for Chemical Solutions brand of PQ Kits ordering information.

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 09-01-2024