



## What is the bending radius of 0.0025 inch ID 1/16 inch OD PEEK tubing - FAQ

The smallest bend for PEEK tubing is a radius of 25 mm which is a diameter of 5 cm (2 inch).

Click [HERE](#) for the ordering information page for PEEK tubing.

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 07-22-2024