



What is the bending radius of 0.0025 inch ID 1/16 inch OD PEEK tubing - FAQ

The smallest bend for PEEK tubing is a radius of 25 mm which is a diameter of 5 cm (2 inch).

Click [HERE](#) for the ordering information page for PEEK tubing.

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com