

What is the best way to track samples in U-2D inserts – Tips and Suggestions

The U-2D™ inserts can be written on using our laboratory markers found [HERE](#).

Also, the U-2D™ racks have positions that are marked with alpha-numeric embossed markings on the sides of the racks.



Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 07-22-2024