



What is the cutting blade material made from on the Clean Cut90 tubing cutter - FAQ

The [Clean Cut90 Tubing Cutter](#) blade is made from standard steel. We wrap this blade in tin foil for safety as well as to prevent rusting when not in use.

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com