



## What is the difference between a Last Drop Filter and the Economy Version - FAQ

The primary difference between the economy version and the Last Drop filter lies in the way the tubing connects to the filter. The economy version features a simple hole for the tubing, which can sometimes result in the filter slipping off. In contrast, the Last Drop filter includes various connectors to securely hold the tubing in place.

Click [HERE](#) for Last Drop Filters™ ordering information and pictures.

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)