



What is the difference between a cap liner and a cap septum - FAQ

Liners or Septa?

A **liner** can be made of the same material and be the same or different thickness as a septum but it is installed in a solid top cap and therefore cannot be pierced as **septa** do. Its primary function is to keep the sample from escaping the vial while the cap is on.

Note: According to Wikipedia a septum (septa is plural) is Latin for something that encloses; or a wall, dividing a [cavity](#) or structure into smaller ones. For chromatography, a septum encloses the vial keeping samples and gases enclosed in the vial until a needle penetrates it and extracts some of the sample to be injected and then reseals itself.

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com