



## What is the temperature range for non-coated Mobile Phase Reservoir Bottles – FAQ

The temperature range for [non-coated Mobile Phase Reservoir Bottles](#) is -70° to 500°C. Please remember to verify the cap's temperature range as well and avoid fully sealing the bottles when heating.

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)