



Does HSQ/PQ software need to be saved as a macro enabled XLSM file - FAQ

YES, the HSQ/PQ software must be macro enabled, the file extension is XLSM.

One small downside of using Excel, lies in its reliance on macro-enabled spreadsheets. Some companies impose stringent restrictions on the use of such macros.

If your laboratory SOP requires the ability to run a spreadsheet, you may need to consult with your IT department to enable this functionality.

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 09-21-2024