

## Wikipedia reference for use of pipets with organic solvents - Primer

### From Wikipedia

#### **Piston-driven air displacement pipettes**

These are the most accurate and precise pipettes. They are more commonly used in biology, though they are commonly used by chemists as well.

The plastic pipette tips are designed for aqueous solutions, and are not recommended for use with organic solvents which may dissolve the plastic.

[Link to Wikipedia page](#)

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)