



## What does the term Activity refer to in chromatography - FAQ

It is defined as the relative strength of the surface of the stationary phase in adsorption chromatography. In the case of silica gel, the more silanol groups are available on the surface, the more active it is. One can modify activity by adding water or other polar solvent that will bond to active sites, reducing activity.

Cogent TYPE-C™ HPLC columns have very low activity due to the virtual absence of silanol functional groups on the surface.



Printed from the Chrom Resource Center  
Copyright 2025, All Rights Apply  
**MicroSolv Technology Corporation**  
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900  
Fax: (910) 769-9435  
Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)  
Website: [www.mtc-usa.com](http://www.mtc-usa.com)