



## The difference between potency and purity in drug products - Tips & Suggestions

**Please refer to the International Conference on Harmonization (ICH) guidelines for more detail.**

A quick definition would be:

-**Potency** is a measure of drug activity expressed in terms of the amount required to produce an effect of given intensity.

-**Purity** is a measure of the amount of API present in a sample compared to those of related substances, impurities, residual solvents, etc.



Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)