


Sample loop information - FAQ

Sample Loops are available in a variety of internal volumes.

Some of these will have several turns of the loop but others will not have a loop at all. In the latter case, the reason is because the tubing length is too small to make a loop. Whether or not there is a loop will not affect the functionality of the product.



Click [HERE](#) for sample loops ordering information and pictures.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com