



What is the volume of Cogent HPLC columns - FAQ

Use this chart for approximate answers to that question or click on the link below for other ID and Length:

<u>Calculate your own</u>	4.6mm ID	Volume
	50mm Long	0.83ml
	75mm Long	1.25ml
	100mmLong	1.66ml
	150mmLong	2.49ml
	250mmLong	4.15ml



Printed from the Chrom Resource Center
Copyright 2025, All Rights Apply
MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900
Fax: (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com