



## Are the Last Drop inlet filters UHPLC compatible - FAQ

**YES.**

MicroSolv's Last Drop™ inlet filters for mobile phase bottles ARE compatible with UHPLC, but it should be noted that for columns with 3µm or smaller particle size packing including UHPLC columns, a 0.2µm pore size filter is recommended. *Flow will be slower due to the smaller pores in the filter.*

Click [HERE](#) for Last Drop Filters™ ordering information and pictures.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)