



## What is best way to know if there is a crack in a CE capillary - FAQ

The best way to know is to look at the entire capillary under either a stereo microscope or using a hand held magnifier offered on our website.

Click [HERE](#) for ordering information and pictures of the magnifier for this purpose.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)