

## What is the difference between a cap liner and a cap septum - FAQ

### Liners or Septa?

A **liner** can be made of the same material and be the same or different thickness as a septum but it is installed in a solid top cap and therefore cannot be pierced as **septa** do. Its primary function is to keep the sample from escaping the vial while the cap is on.

*Note: According to Wikipedia a septum (septa is plural) is Latin for something that encloses; or a wall, dividing a [cavity](#) or structure into smaller ones. For chromatography, a septum encloses the vial keeping samples and gases enclosed in the vial until a needle penetrates it and extracts some of the sample to be injected and then reseals itself.*

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)