



## PEEK column hardware end fitting dimensions - INTERNAL ONLY

- The 7.8 mm PEEK HPLC columns use the end fitting. OD: 20 mm
- 2.1 mm, 4.0 mm and 4.6 mm have an end fitting OD: 12.7 mm (0.5 inches).
- These are standard connection type of ferrules and nuts with 10 / 32 threads.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)