



Does the Mobile Phase Filter made with polyethylene contain PTFE - FAQ

Do the polyethylene Mobile Phase Filters contain PTFE?

These products do not contain PTFE as long as you use a Tefzel / [ETFE ferrule](#), catalog number [58847-01](#). [Mobile Phase Filter, 20 um, Polyethylene 58847-01](#)

[Replacement, Mobile Phase Filter, 20 um, Polyethylene. 58847-02](#)

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com