

What is the best way to track samples in U-2D inserts - Tips and Suggestions

The U-2D™ inserts can be written on using our laboratory markers found [HERE](#).

Also, the U-2D™ racks have positions that are marked with alpha-numeric embossed markings on the sides of the racks.



Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com