



Can I use the MicroSolv Window Maker be used for smaller capillaries - INTERNAL ONLY

INTERNAL

The MicroSolv Window Maker™ is designed for removing a thin film of polyimide from fused quartz silica tubing with 365um OD and ID of 25, 50, 75, 100 and 125um for use in Capillary Electrophoresis. We cannot say if the device will function to your specifications (size 192) with other products as we are not familiar with them.

The Window Maker™ will likely work in this case, as long as one is trying to remove a polyimide coating. The Window Maker could be used with different size capillaries can be used, but the timing of how long to leave in the burner with heat may be different. We do not have a lot experience with smaller capillaries. A definitive answer should not be given.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com