

## Can I backflush the Cogent column filter for longer in-service life - FAQ

### The Cogent column filter can be backflushed but...

Backflushing might push the internal filter element out of place. And not having a clog in the frit, but with slightly dislodged filter element suggests that filter is not working. Therefore, backflushing is not supported.

Click [HERE](#) for ordering information and a picture of the Cogent column filters.



Printed from the Chrom Resource Center  
Copyright 2025, All Rights Apply  
**MicroSolv Technology Corporation**  
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900  
Fax: (910) 769-9435  
Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)  
Website: [www.mtc-usa.com](http://www.mtc-usa.com)