



FlexiTube specifications - Tech Information

The FlexiTube™ length and coil data for selecting the right tube for you.

Ordering information for each type of FlexiTube can be found on the [ordering information](#) page, under each individual item.

Type	Coil resting position	Coil extended position	Coil quantity	FlexiTube length fully extended	Total length of tubing
2	10 mm	30 mm	4	330 mm	710 mm
3	20 mm	120 mm	9	420 mm	945 mm
4	35 mm	200 mm	15	500 mm	1225 mm
5	50 mm	300 mm	23	600 mm	1585 mm
6	100 mm	400 mm	46	700 mm	2060 mm
7	150 mm	500 mm	69	900 mm	3150 mm
8	170 mm	8000 mm	78	1100 mm	4710 mm
9	225 mm	1100 mm	104	1400 mm	6100 mm
10	250 mm	1300 mm	116	1600 mm	6930 mm

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com