



## Physical nature of PEEK frits - Technical Information

Due to the physical nature of the material and how it is manufactured, one must be careful when handling PEEK Frits.

PEEK frits can be brittle. Do not throw the bag of frits or drop it. Use caution when removing the frits in the bag, as bumping into one another can cause dislodging of the material.

Printed from the Chrom Resource Center  
Copyright 2025, All Rights Apply  
**MicroSolv Technology Corporation**  
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900  
Fax: (910) 769-9435  
Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)  
Website: [www.mtc-usa.com](http://www.mtc-usa.com)