



Ferrule selection suggestions - Tech Information

When using a fitting nut, tubing and ferrule, the ferrule should be the softest material so that it can swage on to the tubing and hold it in place and be "formed" by the nut.

But always use PEEK ferrules with PEEK tubing.

[More information about PEEK tubing](#)

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com