



Which ferrules are best when using soft polymeric tubing - FAQ

Kel-F™ (*PCTFE which is similar to PTFE*) and Tefzel™ (*ETFE*) ferrules give the best results when using PTFE, FEP, PFA or ETFE tubing.

[More information about PEEK tubing](#)

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com