



What can I use for a reference solution against Holmium Oxide - FAQ

Water is an appropriate reference solution. Since the lowest wavelength in Wavelength Calibration Solution is 241 nm, any substance that doesn't absorb in that region can be used as a reference solution.

Click [HERE](#) for Chemical Solutions brand of PQ Kits ordering information.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com