

Changing the filter element in the Last Drop Filters. - FAQ

We do not recommend changing the filter yourself. During the process, the part holding the filter element can be damaged, potentially creating a bypass or causing the frit to no longer stay in place.

Click [HERE](#) for Last Drop Filters™ ordering information and pictures.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com