



## What is the temperature range for non-coated Mobile Phase Reservoir Bottles - FAQ

The temperature range for [non-coated Mobile Phase Reservoir Bottles](#) is -70° to 500°C. Please remember to verify the cap's temperature range as well and avoid fully sealing the bottles when heating.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)