



## Is the Zero and Controlled flow capillary available in larger than 1-meter length - FAQ

Due to the difficulty of properly bonding the inner walls of the Zero and Controlled flow capillaries, they are only available in 1-meter lengths.



Click [HERE](#) for ordering information and pictures

Printed from the Chrom Resource Center  
Copyright 2025, All Rights Apply  
**MicroSolv Technology Corporation**  
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900  
Fax: (910) 769-9435  
Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)  
Website: [www.mtc-usa.com](http://www.mtc-usa.com)