



Can the all silicone cap mats be sterilized - FAQ

Yes, the silicone cap mats can be sterilized using heat at 121°C for 15-30 minutes at 16-17 psi.

Click [HERE](#) for Silicone rubber sealing mat ordering information and pictures.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com