



What are the temperature specifications of the autosampler vial storage box and cap for 4ml vials - FAQ

Catalog number 9504SB-03 is temperature resistant in the range -80°C to +100°C.

Click [HERE](#) for storage box for 15x45mm autosampler vial ordering information.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com