



What is the internal diameter of the opening for 8mm crimp top BASIK brand autosampler vials - FAQ

The ID of the opening for 8mm crimp top MicroSolv BASIK vials, catalog numbers [9508C-FB08-CV](#), [9508C-FB12-CV](#), and [9508C-FB12-AV](#) is **4.0mm +/-0.2**.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com