



## What is the advantage of using a Pharmafix type of septa for headspace analysis - FAQ

**The advantage is that you get a better seal with the “Pharmafix” seal.** However, you should be aware that it may not be as inert as a silicone/PTFE septa designed for headspace analysis.

Click [HERE](#) for ordering information and pictures of headspace caps

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)