



Which crimp cap is most commonly used for residual solvent testing by the pharmaceutical industry - FAQ

The crimp cap with silicone rubber / PTFE septa is most used by the pharmaceutical industry to be in compliance with USP <467>.

Our catalog number [95025-09-1S](#) is one of many caps used. Click [HERE](#) for more information about Headspace vials and caps

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com