



For the bonded or coated CE capillaries do you have other lengths available - FAQ

The bonded or coated capillaries are only available in 1M lengths.

The reason is that it is very difficult to push the bonded reagent through a 50um or 75um tube and get proper coating. Click [HERE](#) for ordering information Controlled Flow and Zero EOF Capillaries

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com