



Do you sell CE capillaries of various outer diameters - FAQ

No, the only outer diameter (OD) bare fused silica capillary for CE that we sell is 365µm.

Click [HERE](#) for capillaries for CE ordering information and pictures.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com