



Are there other parts or pieces that I need to use the MicroSolv Window Maker - FAQ

No, the MicroSolv Window Maker™ comes with a heating module that makes a 2mm window. No other tools are needed. Other heating modules of different window sizes can be purchased separately.

Click [HERE](#) for MicroSolv Window Maker ordering information and pictures.

Printed from the Chrom Resource Center
Copyright 2025, All Rights Apply
MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900
Fax: (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com