

What is the polyimide cladding on CE capillaries - FAQ

Polyimide is an aromatic, linear polymer.

It is not a true thermoplastic because it thermally degrades before its glass transition temperature. It is not cross linked. It is used on CE capillaries because it is solvent resistant, possess barrier properties and performance at both high and low temperatures.

Click [HERE](#) for Capillary Technical Properties.

Click [HERE](#) for MicroSolvCE Window Maker ordering information and pictures

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com