


## How do I stack multiple ColumnStore column storage cabinet units - FAQ

When ordering multiple units of the The ColumnStore™ column storage boxes are shipped with labels placed on each ColumnStore unit which read "Box 1 of 3"/ "Box 2 of 3"/ "Box 3 of 3" for example. This will indicate the order that the boxes should be stacked.

Follow the instructions in the below to stack the boxes. You must use a mounting bracket (special order item) to join two units together. The top box does not have a mounting bracket.

If you order the units at different times, simply order the mounting bracket shown below and stack the units placing the bracket between the storage units.



Printed from the Chrom Resource Center  
Copyright 2025, All Rights Apply  
**MicroSolv Technology Corporation**  
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900  
Fax: (910) 769-9435  
Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)  
Website: [www.mtc-usa.com](http://www.mtc-usa.com)